

The RX7000 Alpha Automatic Digital Refractometer



Features:

- Brix Range 0.0 – 95%
- Refractive Index 1.3250 – 1.7000
- Internal Temperature Control 5 - 70°C (no water bath required)
- Flash memory for 30 values
- Up to 30 user scales
- RS232 Interface



Specifications:

Measurement range.....	Brix 0.0 to 95.0%
Resolution.....	0.01% Brix, 0.00001 RI, 0.01°C
Measurement Accuracy.....	+/- 0.1% Brix, +/- 0.00001 RI, +/- 0.05°C
Measurement Temperature.....	5.0 to 70°C (Automatic temperature compensation)
Prism.....	Artificial Sapphire
Dimensions and Weight.....	37 x 26 x 14 cm, 9.0 kg

Applications:

Alcoholic Drinks	Puree	Animal Oils and Fats
Fruit Drinks	Canned Syrup	Edible Oils
Soft Drinks	Condensed Milk	Vegetable Oils and Fats
Tea Drinks	Ketchup	Flavors
Fruit concentrates	Soup	



WILKENS-ANDERSON COMPANY

4525 West Division Street, Chicago, Illinois 60651
Phone:773-384-4433, 800-847-2222, FAX:773-384-6260