



WILKENS-ANDERSON COMPANY

AUGUST SPECIAL

Special sale on all disposable transfer pipettes through 9.30.15!

Samco Disposable Transfer Pipettes

- Grips easily with latex gloves - no slipping
- Clear graduated markings
- Fine tip accuracy
- Packaging materials greatly reduce airborne mold spores from entering testing areas
- Inert, nontoxic Polyethylene
- Safe, residue-free product that does not interfere with testing data
- Expert craftsmanship produces largest variety of fine tip pipettes in the industry
- Precision-cut tips deliver a uniform drop size every time
- Proprietary controlled processes produce a uniform bulb with uncompromised integrity

Product Specifications

Key #	WACO Product Number	Samco Number	Item Description	Length	Overall Capacity	Bulb Draw	Drops per mL	Qty per Pack	Qty per Case	Samco List Price per Pack	WACO Sale Price per Pack
General Purpose Transfer Pipettes											
1	71603-01	202	General Purpose Standard	15.0 cm	7.0 mL	3.2 mL	22	500	5000	\$30.70	\$23.50
2	71603-07	204	General Purpose Large Bulb	15.0 cm	8.0 mL	4.6 mL	22	400	4000	\$28.30	\$21.68
3	71603-08	273	General Purpose 3"	8.7 cm	1.7 mL	0.9 mL	25	500	4000	\$30.00	\$24.73
4	71603-09	335	5" B/B-Pet (Blood Bank)	12.9 cm	4.0 mL	1.9 mL	25	400	4000	\$29.00	\$22.21
5	71603-05	336	6" B/B-Pet (Blood Bank)	15.5 cm	5.0 mL	1.9 mL	25	500	5000	\$28.30	\$24.51
Graduated Transfer Pipettes											
7	71603-02	212	Graduated 1 mL Small Bulb	13.8 cm	3.0 mL	2.1 mL	25	500	5000	\$28.70	\$21.97
8	71603-03	222	Graduated 1 mL Large Bulb	15.0 cm	5.0 mL	3.1 mL	25	500	5000	\$31.30	\$23.97
8	71603-13	222-20S	Packs of 20, Sterile	15.0 cm	5.0 mL	3.1 mL	25	500	5000	\$51.30	\$36.23
8	71603-93	222-1S	Ind. wrapped, Sterile	15.0 cm	5.0 mL	3.1 mL	25	500	5000	\$70.40	\$49.70
9	71603-04	225	Graduated 3 mL Large Bulb	15.0 cm	7.0 mL	3.2 mL	22	500	5000	\$28.30	\$21.68
10	71603-06	274	Graduated 2 mL	15.6 cm	5.0 mL	1.8 mL	25	500	5000	\$31.70	\$22.38
Fine Tip Transfer Pipettes											
12	71603-30	232-11	Fine Tip Standard	14.4 cm	5.0 mL	3.3 mL	50	500	5000	\$53.00	\$37.43
13	71603-31	234	Fine Tip Large Bulb	15.0 cm	8.0 mL	4.3 mL	45	400	4000	\$43.30	\$30.57
14	71603-35	231	Extended Fine Tip Small Bulb	10.4 cm	1.5 mL	1.0 mL	50	N/A	3200	\$350.40	\$212.59
Narrow Stem Transfer Pipettes											
21	71603-20	241	SEDI-PET Short Stem	6.3 cm	1.2 mL	0.9 mL	25	500	4000	\$31.60	\$22.32
23	71603-15	251	SEDI-PET Standard	15.0 cm	4.0 mL	3.1 mL	25	500	5000	\$33.00	\$25.27
25	71603-17	694	Q-Pet Transfer Pipet	15.5 cm	6.5 mL	3.0 mL	N/A	500	5000	\$40.70	\$28.74
Extra Long Transfer Pipettes											
26	71603-18	262	Extra Long 9"	22.5 cm	6.0 mL	2.3 mL	25	400	4000	\$31.70	\$24.27
27	71603-10	263	Extra Long 12"	30.0 cm	23.0 mL	7.1 mL	22	100	1000	\$36.00	\$25.42