



Bellingham + Stanley Ltd.

**The RFM 100/300 Series
Automatic Digital Refractometers**



- Artificial sapphire prism
- Easy wipe stainless steel sample dish
- External power supply
- RS232C interface
- Large LCD display
- User friendly software
- Password security (two level)
- Optional printing

Specifications	RFM110	RFM130	RFM311	RFM320	RFM330	RFM340
Scales	0-50% Brix	0-95% Brix	0-50% Brix 1.33-1.42 RI User Defined	0-50% Brix 1.33-1.42 RI User Defined	0-95% Brix 1.33-1.53 RI User Defined	0-95% Brix 1.33-1.53 RI User Defined
Resolution	0.1% Brix	0.1% Brix	0.1% Brix 0.0001 RI	0.01% Brix 0.00001 RI	0.1% Brix 0.0001 RI	0.01% Brix 0.00001 RI
Accuracy	+/- 0.1% Brix	+/- 0.1% Brix	+/- 0.1% Brix +/- 0.0001 RI	+/- 0.03% Brix +/- 0.00004 RI	+/- 0.1% Brix +/- 0.0001 RI	+/- 0.03% Brix +/- 0.00004 RI
Reading Time	3.5 sec	6.5 sec	3.5 sec	3.5 sec	6.5 sec	6.5 sec
List Price	\$4,525.00	\$6,180.00	\$7,540.00	\$8,344.00	\$8,089.00	\$9,648.00

The RFM Series of Digital Automatic Refractometers from Bellingham + Stanley meets the requirements of highest precision and repeatability for a wide range of industrial applications of refractometry. All instruments in the range derive from highly developed designs utilising state-of-the-art technologies and are crafted to the highest quality standards with the best available materials. They have achieved a proven track record with customers around the world and have been shown to be ideally suited to many quality control operations.



WILKENS-ANDERSON COMPANY

4525 West Division Street, Chicago, Illinois 60651

Phone: 773-384-4433, 800-847-2222, FAX: 773-384-6260