



Bellingham + Stanley Ltd.



**Introduces the
NEW RFM 700 Series**

- Accurate
- Reliable
- Robust
- Durable
- Moisture Resistant
- Simple to Use
- RS-232
- Low cost

Specifications	RFM710	RFM730
Scale	% Sugar (° Brix)	
Range	0 – 50	0 – 95
Resolution	0.1	0.1
Accuracy	+/- 0.1	+/- 0.1
Reading Time	3.5 seconds	6.5 seconds
Working Temperature Range	5 – 40 °C	5 – 40 °C
Automatic Temperature Compensation	To 20 °C	To 20 °C

RFM700 refractometers are low cost, fully automatic instruments that are equipped with the very latest optical and electronic technology from Bellingham + Stanley. They are ideally suited to quality control in food and beverage industries, but can also be used in many other applications. The no-moving-parts design ensures reliable and repeatable readings. RFM700 instruments are designed for simple operation via the four keys:

- **Read** – gives direct reading in Brix; also power standby function
- **Zero** – allows zero point calibration with water
- **Temperature** – displays temperature in °C
- **Print** – print results

Call us for information and prices

WACO

WILKENS-ANDERSON COMPANY

4525 West Division Street, Chicago, Illinois 60651

Phone:773-384-4433, 800-847-2222, FAX:773-384-6260