



Bellingham + Stanley Ltd.

**Introduces the
NEW RFM 800 Series**

Internal temperature control
(No water bath required)

- Intelligent Thermal Management
- Rapid sample stabilization
- Sample temperature control up to 100 degrees
- New optical system
- High accuracy



Specifications	RFM830	RFM840	RFM870
Scales RI Sugar (%Brix)	1.32-1.58 0-100	1.32-1.58 0-100	1.30-1.70 0-100
Resolution RI Sugar (%Brix)	0.0001 0.1	0.00001 0.01	0.00001 0.01
Accuracy RI Sugar (%Brix)	+/- 0.0001 +/- 0.1	+/- 0.00002 +/- 0.02	+/- 0.00002 +/- 0.02
Automatic Temperature Compensation	5 – 80 °C		
Reading time (s)	Adjustable for accuracy (3-10 sec)		
Temperature Sensor Accuracy	+/- 0.03 °C		
Prism	Artificial Sapphire		
Prism Dish	316 Stainless Steel		
Temperature Control	Peltier		

Call us for information and prices

WACO

WILKENS-ANDERSON COMPANY

4525 West Division Street, Chicago, Illinois 60651

Phone:773-384-4433, 800-847-2222, FAX:773-384-6260