

The Reichert AR200 Digital Hand-Held Refractometer

Full Range Brix and Unparalleled Accuracy in the Palm of Your Hand

The Reichert logo features a stylized blue circle with a white shape inside, resembling a drop or a lens, positioned above the word "Reichert" in a bold, blue, sans-serif font.

- Digitally measure refractive index from 1.3330 to 1.5600 or % solids as Brix over the full 0 – 95% range
- Measures transparent, translucent and opaque samples
- Achieve exacting accuracy to +/- 0.1% Brix or 0.0001nD with 3600 cell, linear scanned array imaging
- Exclusive IR communications port allows easy downloading of readings directly into your PC to meet compliance needs (optional accessories required)
- Prism temperature sensor and automatic temperature compensation enhance accuracy without temperature correction calculations (10-110°F/-12-43°C)
- 100's of custom channels for use in beverage, food, chemical and agricultural industries are available or can be created to meet your analytical needs.



List Price.....\$1099.00

**WACO Special
Price.....\$995.00**

The WACO logo consists of the letters "WACO" in a large, bold, red-outlined font.

WILKENS-ANDERSON COMPANY

4525 West Division Street, Chicago, Illinois 60651

Phone:773-384-4433, 800-847-2222, FAX:773-384-6260