

# TempTestr IR Noncontact Temperature Measurement

**OAKTON**  
I N S T R U M E N T S

- **Infrared measurement provides safe, noncontact readings** – measure moving targets or targets too hot to touch; avoid contamination problems
- **Simple one-button operation** – just point and press
- **Laser sighting** – pinpoint your exact target area
- **Extremely fast results** – only ½ second for 95% response
- **Hold function** – great for hard to reach areas
- **Selectable °F/°C display**
- **Ergonomic, compact design** – for easy one handed operation
- **-25 to 400°F / -30 to 200°C range**
- **+/-1°C (2°F) accuracy** – at HACCP critical temperature range of 0 to 60°C/32 to 150°F
- **Meets UL standards for IR thermometers in the food industry**



**WACO**

---

## **WILKENS-ANDERSON COMPANY**

4525 West Division Street, Chicago, Illinois 60651  
Phone:773-384-4433, 800-847-2222, FAX:773-384-6260

---